

Presumptive Drug Identification Using Chemical Reagents

1. Purpose and scope: This procedure outlines standardized steps for conducting **presumptive drug identification** using chemical reagents in field conditions.

- **Audience:** Border Guard Forensic Technicians.
- **Environment:** Roadside inspections and border crossings.
- **Constraint:** Results are **presumptive only** and require accredited laboratory confirmation.



- **Safety precautions:**
 - **Wear:** gloves, goggles, and a protective apron
 - **Ventilation:** Conduct tests in well-ventilated areas
 - **Handling:** Handle reagents with care; avoid contact with skin/eyes
 - **Disposal:** Follow hazardous waste protocols.

3. Required Materials

- **Glassware:** Transparent vials or test tubes with stoppers.
- **Instruments:** Glass droppers/pipettes and a clean glass rod for crushing.
- **Kits:** Field reagent kit and waste container

Instructions

4. Sample Preparation

- Place a pinch-sized quantity of the substance into a clean, dry vial.
- For tablets or plant material: Gently crush using a clean glass rod.

 **Important:** Do **not** add any neutral liquids; the reagent is the sole reactive medium.

5. Reagent Reference Table

Reagent	Target Substances	Expected color change
Marquis	Opiates, amphetamines	Purple (opiates), orange-brown (amphetamines)
Mandelin	Methamphetamine, ketamine, MDMA	Green (meth), black (MDMA)
Mecke	Heroin, morphine	Blue-green to black
Scott Test	Cocaine	3-step: Blue → Pink → Blue (layered)
Duquenois–Levine	Cannabis	Purple in the chloroform layer
Simon's	Secondary amines (e.g., MDMA)	Blue

6. Testing Procedure

- **Application:** Add **1–2 drops** of reagent directly onto the sample.
- **Selection:** Choose the reagent based on the suspected substance.
- **Observation:** Monitor for up to **60 seconds** until a color change occurs.

7. Interpretation and Documentation

Compare results to the kit's reference chart and record the following:

- Tested substance and reagent used.
- Time elapsed to change and final color observed.
- Reagent batch number

 **Final Note:** Mark results clearly as "**presumptive only**".

8. Maintenance and Disposal

- **Neutralization:** Use baking soda solution before disposal.
- **Cleaning:** Rinse glassware with distilled water and air dry.
- **Storage:** Store reagents in sealed, labeled containers away from heat/light.

 Powiązane artykuły [Forensic Conference at BGTC](#)

Filtruj wg etykiety

Teraz nie ma elementów z wybranymi etykietami.